



ActivThumb™ is easy to wear and can be used following fractures, injuries to the thumb and for skiers thumb.




ActivThumb™ supports the MP joint (where the thumb joins the hand) and aids common problems such as arthritis (osteoarthritis) and instability. Decreases pain by stabilising the thumb and helps to prevent any further injury.

Built in cross belts are used to apply strong fixing. Cushioning material and anti-slip allow this product to be worn during sporting activities as other functions of the hand are not limited.


How to wear:



Loosen the cross belt and open the thumb hole



Insert your thumb and close the secondary velcro fastener



Pull the cross belts tight and secure them



Secure the main velcro fastener

Caution

- When wearing for the first time, try this product for 1 to 2 hours
- Do not wear while sleeping
- Do not use this product other than explained in these user instructions as this can lead to damages or injuries
- make sure the fastener does not stick to clothes as this may cause damage to clothes
- Do not close the belt too tightly. Blood circulation and skin damage may occur if product is used for prolonged periods
- This product contains chemical fibres. If the product irritates the skin stop wearing the product immediately and consult your doctor

Washing:

- Hand wash with cold water
- No whitening or bleaching ingredients
- Air dry with good ventilation away from direct sunlight
- This product may lose colour because of water, sweat or friction

Storage:

- Store in a dry place with good ventilation
- Do not expose to direct sunlight or high temperature

Disposal:

- Please conform with local guidelines for textile product disposal

ActivThumb™ has a removable flexible plastic stay and cross belt mechanism to ensure a good fit and give you the confidence you need to continue your activities. This product is perfect for sports people who need to continue using their hand e.g. basketball and volleyball players.





Colour: Black

Circumference: 13 - 21 cm

Unisex

Handmade in Japan by:

 Daiya Industry Company, Ltd.
 1117-1 Koshinden, Minami-Ku Okayama 701-0203 Japan

Distributed in UK, EU and Australia by:

 **Gower**
Health & Fitness
Solutions

271 High Street, Berkhamsted HP4 1AA
United Kingdom
info@activsupports.com